

EASIMOVEPRO™

EVERY **TRANSFER** MATTERS.REUSABLE AIR-ASSISTED
LATERAL TRANSFER MATTRESSES

EasiMove™ PRO supports reusable lateral transfer workflows for facilities looking to reduce reliance on single-use transfer products while maintaining familiar air-assisted patient handling practice.



01

Reusable Value

Long-term patient handling value



Reduced single-use reliance. Supports facilities seeking to reduce ongoing consumable use and associated clinical waste.



Reusable workflow model. Designed for repeated clinical use with appropriate cleaning, maintenance and inspection.



02

Heat-Sealed Construction

Cleanability and repeated-use durability



Heat-sealed seams. Helps reduce stitch-line and needle-hole concerns associated with sewn construction.



Wipe-clean, launderable design. Supports routine wipe-down between uses and machine laundering when required.



03

Clinical Simplicity

Set-up, alignment and compatibility




Visual positioning cues. Head position indicator and centre line support patient orientation, centring and set-up.



Air-supply compatibility. Compatible with EasiAir™ Variable Speed Air Supply and HoverTech® HT-Air® 2300 Air Supply.

PRODUCT OVERVIEW

EasiMove™ PRO is available in 34" and 39" widths to support reusable air-assisted lateral transfer workflows across compatible care surfaces.



KEY 1. Printed Head Outline 2. Patient Safety Straps 3. Reinforced Transfer Handles 4. Centre Alignment Line
5. Product Identification Label 6. Inflation Port Connections

SIZES & SPECIFICATIONS

SPECIFICATION	PRO 34"	PRO 39"
Product code / REF	MHDM-EM34PRO	MHDM-EM39PRO
Width	86 cm	99 cm
Length	198 cm	198 cm
Construction	RF-welded / heat-sealed	RF-welded / heat-sealed
Material	Nylon twill	Nylon twill
Reusable	Yes	Yes
Safe working load	453 kg	453 kg
Air supply compatibility	EasiAir™ / HT-Air® 2300	EasiAir™ / HT-Air® 2300

SUMMARY

The EasiMove™ PRO range combines reusable workflow value, heat-sealed construction, point-of-use guidance and air-supply compatibility within one familiar reusable lateral transfer platform.

Every **Move** Matters.

HoverTech® and HT-Air® are registered trademarks of HoverTech International. Product names are used for comparative reference only.